



PRODIGY NPWT Device



“The Intelligence to Heal”

Main Features:

- **Controlled by micro processor**
- **Touch screen operation**
- **Maintenance free device**
- **Adjustable vacuum pressure**
- **Altitude compensation**
- **AC power with battery backup**
- **Alarms for:**
 - Canister Full
 - Blocked Filter
 - Air Leak
 - High Vacuum
 - Low Battery

699 Main Street, New Rochelle, NY 10801
(914) 632-3988 Fax: (914) 636-8852
www.premcomedical.com